

# ACCU-Beam®

# TTI Medical

## *ACCU-Beam® Digital SLR Camera Adaptor for Microscopes and Slit Lamps....*



### *FEATURES*

- **For Canon, Nikon, Fuji, Sigma, Kodak, Olympus, Pentax and Minolta SLR Digital Cameras**
- **Attaches directly to Zeiss® style Beamsplitters**
- **Beamsplitter adaptors available for Leica, Global Moller-Wedel, Kaps, Topcon and Haag-Streit 900BQ**
- **Convenient Fine Focus Control**
- **Precise Image Centering Adjustment**
- **Adjustable Iris to control Depth of Focus and Image Brightness**

# ACCU-Beam® Digital SLR Camera Adaptor

## Specifications:

Optics	AR coated
Transmission	≥90%
Focal Length	F200
Focus Control	+/- 2mm
Iris Control	adjustable 0-15mm
Centering	X-Y - full range
Dimensions	5.5" L x 4.5" H
Weight	20 oz

**8194101 Canon EOS**

**8194102 Nikon, Fuji, Sigma, Kodak**

**8194103 Minolta**

**8194104 Olympus**

**8194105 Pentax**



## Canon EOS 6.3 MP Digital Rebel with Custom Footswitch

### Beamsplitter Adaptors:

8110K Kaps	8115 Topcon SL2E, 4E, 6E, 2D, 5D
8111 Global	8116 Topcon OMS-300, 310, OMS-60, 70
8112 Leica / Wild M650, M690	8117 Haag-Streit 900-BQ
8112/R Leica OHS, MS, M500	8118 Rodenstock
8113/2 Moller-Wedel	

TRANSAMERICAN TECHNOLOGIES INTERNATIONAL

2246 Camino Ramon • San Ramon, CA 94583 • USA

Tel: 925-355-0750 • 800-322-7373 • Fax: 925-355-0777

e-mail: [info@ttimedical.com](mailto:info@ttimedical.com) • website: [www.ttimedical.com](http://www.ttimedical.com)



ISO 9001 / EN 46001 / ISO 13485



Marked Products