


ACCU-Beam®

TTI Medical

***ACCU-Beam® Universal CO₂ Laser
Handpieces...***




ACCU-Beam® CO₂ Laser Handpieces are available in focal lengths of 75mm, 100mm, 125mm and 200mm. In addition there is a variable spot handpiece for Dermatology and a zoom handpiece for GYN/Procto procedures.

Each handpiece is universally adaptable to all commercial articulated arm CO₂ lasers by using articulated arm thread adaptors. Adaptors can be customized to fit any thread specification on request.

The lens assembly is removable allowing the handpiece to be sterilized by steam autoclave, ETO or cold soak sterilization.

ACCU-Beam® Universal CO₂ Laser Handpieces



6150 GYN/PROCTO ZOOM F250-450

CO₂ Laser Articulated Arm Thread Adaptors

1101	Sharplan Quick Disconnect for 1040, 1060 and 1100
1102	LaserSonics / NIIC
1103	LaserSonics/Merrimack
1104	Coherent / Xanar
1106	Coherent 450 and 451
1108	Zeiss
1109	Sharplan Twist-Lock for 1020, 1050, 1055, 1075 and Ultra Pulse
1112	LaserSonics LS-500

Note: Sharplan 1060 has two arm versions – 1101 & 1109. Thread adaptors are not required for Surgilase, laser Engineering, Sharplan 720, 733A, 734 and LaserSonics Illumina 40 Black Arm.

ISO9001:1994
AMTAC Cert. No. 163



TRANSAMERICAN TECHNOLOGIES INTERNATIONAL

2246 Camino Ramon • San Ramon CA 94583 USA

Tel: 925-355-0750 • 800-322-7373 • Fax: 925-355-0777

e-mail: ttmed@aol.com • website: www.ttimedical.com

Rev:060399