


## *ACCU-Beam® Universal Beamsplitter EndoAdaptors for Rigid and Flexible Endoscopes....*



### *Features*

- **Focal lengths - 25mm and 30mm**
- **Convenient direct viewing eyepiece**
- **Balanced on axis camera mount**
- **Finger tip focus control**
- **Free spinning eyepiece clamp with lock**

# ACCU-Beam® Universal Video Beamsplitter EndoAdaptor

8005 25mm Beamsplitter EndoAdaptor  
8006 30mm Beamsplitter EndoAdaptor



## ***Specifications:***

Focus	Adjustable focus ring
Camera Mount	C-mount
Eyepiece Clamp	Fits 32mm rigid scope eyepieces lockable and free spinning
Soak Disinfecting	Use 2% activated glutaraldehyde
Sterilization	ETO, Steris or Sterrad
Dimensions	1.5 x 2.0 inches
Weight	4.5 oz

**Flexible Endoscope Adaptors:** Connects flexible scope eyepiece to eyepiece clamp.

8050 Pentax  
8051 Machida  
8052 Olympus OES



ISO9001:1994  
AMTAC Cert.No.163



TRANSAMERICAN TECHNOLOGIES INTERNATIONAL  
2246 Camino Ramon • San Ramon, CA 94583 • USA  
Tel: 925-355-0750 • 800-322-7373 • Fax: 925-355-0777  
e-mail: [info@tmedical.com](mailto:info@tmedical.com) • website: [www.tmedical.com](http://www.tmedical.com)