


ACCU-Beam®

TTI Medical

ACCU-Beam® Universal Video EndoAdaptor for Rigid and Flexible Endoscopes....



The ACCU-Beam® Universal Video EndoAdaptor is available in focal lengths of 24mm, 30mm and 35mm. Other focal lengths can be customized on request.

A C-mount thread is used for attaching standard C-mount video cameras. The EndoAdaptor can be modified to accept cameras with other thread mounts.

The eyepiece clamp accommodates a standard 32mm eyepiece. The eyepiece can be locked in place or unlocked to rotate the scope within the clamp. Eyepiece adaptors are available for a wide variety of flexible endoscopes.

ACCU-Beam® Universal Video EndoAdaptor

8009	24mm EndoAdaptor
8012	30mm EndoAdaptor
8013	35mm EndoAdaptor

Specifications:

Focus	Adjustable focus ring
Camera Mount	C-mount and custom camera mounts
Eyepiece Clamp	For 32mm rigid scope eyepieces – lockable
Soak Disinfecting	Use 2% activated glutaraldehyde
Sterilization	ETO, Steris or Sterrad
Length	1.75 inches
Weight	2.4 oz

Flexible Endoscope Adaptors: Connects flexible scope eyepiece to eyepiece clamp.



8050	Pentax	8051	Machida
8052	Olympus OES	8053	Olympus pre-OES
8054	Fujinon	8055	Schott LS7
8056	Schott FS100 & FS101		

ISO9001:1994
AMTAC Cert.No.163



TRANSAMERICAN TECHNOLOGIES INTERNATIONAL
2246 Camino Ramon • San Ramon, CA 94583 • USA
Tel: 925-355-0750 • 800-322-7373 • Fax: 925-355-0777
e-mail: info@ttimedical.com • website: www.ttimedical.com